

## **DNA SAMPLE DELIVERY INSTRUCTIONS**

### **Blood sample's extraction and storage conditions**

- Blood sample of 10 ml.
- Complete blood count in EDTA tube.
- No centrifugation needed.
- Store at room temperature. Never freeze.

### **Sending the blood sample**

The blood sample must be sent through a transport company to the following address:

**Felix Claverie-Martin**  
**FEA Unidad de Investigación (2ª Planta de Traumatología)**  
**RenalTube**  
**Hospital Universitario Nuestra Señora de Candelaria**  
**Carretera del Rosario 145**  
**38010 Santa Cruz de Tenerife**  
**España**

Proper actions must be taken so that the blood sample arrives between 9:00h - 14:00h Monday to Friday.

Blood samples must be accompanied by patient's basic data and the referring hospital official document assuming the conformity with the genetic study procedure cost.

Before sending the sample, it is imperative to contact the Research Unit of Hospital Nuestra Señora de Candelaria:

**Félix Claverie Martín**

**Phone: (+34) 922 60 05 66**

**Correo electrónico: [fclamar@gobiernodecanarias.org](mailto:fclamar@gobiernodecanarias.org)**